1. Post-Treatment Preparation
☐ Begin instrument processing immediately after patient care
☐ Keep contaminated instruments contained and secure
 2. Personal Protective Equipment (PPE)
□ Wear nitrile gloves
□ Wear protective eyewear
☐ Wear gown and mask during transport
3. Pre-Cleaning
☐ Ensure instruments and bur blocks are visibly clean
□ Remove all bioburden (blood, tissue, debris)
☐ Use ultrasonic cleaner for thorough cleaning
☐ Do not place dirty instruments directly into sterilization pouches
4. Inspection Before Packaging
□ Inspect instruments for residual debris
□ Check for damage or wear
0
□ Only package clean, intact instruments for sterilization

Pre-Sterilization Protocol Checklist